

PORTABLE CONTACT ANGLE ANALYSER PGX

The PGX is a compact instrument that measures wetting and absorption on virtually any size or shape surface without the need to cut sample strips required by other contact angle testers.

This enables the surface testing of almost any 3dimensional product including :

- a table top
- the roof or bumper of a car
- drywall board
- a glass bottle
- a metal can

The PGX provides automated testing.

Simply place the instrument on top of your specimen and press a button. The internal pump delivers a precise droplet, which is automatically delivered on the test surface. The built-in camera captures a video sequence of the liquid droplet and results are presented as static contact angle or dynamic wetting and liquid penetration as a function of time. These results are invaluable for evaluating and controlling a wide range of processes including adhesion, corona treatment, absorption, surface contamination, surface sizing, printability and wettability.

PGX computer control

The PGX offers the ease of computer-controlled testing, analysis and reporting. Simply connect the PGX to a standard laptop or PC with an USB interface and start the software.

The PGX receives power from the USB connection and a live video image will appear in the computer monitor.



The PGX provides contact angle of any size specimen

Features

- Measurement of :
- Static contact angles wetting of non-sorptive surfaces
- **Dynamic contact angles** wetting, absorption and spreading as a function of time
- Surface tension & surface free energy evaluates the probing liquid from a pendant drop
- · Lightweight, portable unit
- Small foot print : 55 x 90 mm
- · Simple, accurate calibration routine
- Automatic droplet application
- Integrated camera : captures 15 images/sec
- Built-in pump : adjustable in 0.5 µl steps
- Easy to install no hardware ! Runs on most PC's and laptops (USB)
- User-friendly software for Windows 98/ME/2000/XP (with dialogue in your own language)

Let the PGX measuring head X-ceed your X-pectations

Physical specifications

Dimensions 9 x 5.5 x 9 cm (LxWxH)

Net Weight 400 grams

Standards


TAPPI T-458, ASTM D-724, D-5946

Applications

Printing - Agriculture - Paper - Detergents - Corona treatment - Coating - Pharmacy - Board - Inks - Flame treatment - Adhesion - Biomedicine - Wood - Oils -Surface sizing - Cleaning - Papermaking - Metals -Surfactants - Surface tension - Absorption - Cosmetics - Plastic - Water -Surface contamination - Wettability -Dentistry - Textile - Solvents - Surface energy




Nijverheidslaan 47 8540 Deerlijk Belgium Tel.: +32(0)56.78.21.70 Fax: +32(0)56.77.30.40 Email: info@rycobel.be Website: www.rycobel.be




rycobel*group*







1. Static results from 5 + 4 droplets on two solid substrates





2. Dynamic results from a droplet on an absorbent substrate





3. Stored data can be retrieved by the Viewer or a spreadsheet




















iquid: Undefined	Density: 1.0				Measu
	-1				द्भुभू Measu
85	÷				
60					
55					
50					
45					
40					
30					
25					
20					
15					
5					
°.[
0	1 2	3 4	5 6	Avg 67.71	New
	Surface Tension				
T	71.13 67.52				🔂 Oper
	67.34				
	67.04 67.77				🔒 Save
	66.79				
	66.38				🕕 Exit
Dynamic Advancing	Pol Dochark Till Liferenenf Frogram version 2-2 SuperLow View Surface Te	nsion Surface Energy Calibrate H	elp		
Dynamic Advancing	SuperLow View Surface Te		elp		
and the second	SuperLow View Surface Te	🔗 Browse	elp		
Dynamic Advancing uple ID:	SuperLow View Surface Te		elp		
Dynamic Advancing ple ID: Culation Model ASTM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water:	🔗 Browse	eip		
Dynamic Advancing ple ID: sulation Model STM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse	elp		-
Dynamic Advancing nple ID: culation Model ASTM D5946	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse	elp		
Dynamic Advancing ple ID: Culation Model ASTM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse	elp		
Dynamic Advancing ple ID: culation Model ASTM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse	eip		
Dynamic Advancing ple ID: culation Model ASTM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse	elp		-
Dynamic Advancing ple ID: culation Model ASTM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse	eip		
Dynamic Advancing ple ID: Culation Model ASTM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse	elp		
Dynamic Advancing ple ID: Culation Model ASTM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse	eip		
Dynamic Advancing ple ID: Culation Model ASTM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse	elp		
Dynamic Advancing ple ID: culation Model ASTM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse	eip		
Dynamic Advancing ple ID: Culation Model ASTM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse			Calc
Dynamic Advancing ple ID: culation Model ASTM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse	elp		Calc
Dynamic Advancing ple ID: culation Model ASTM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse			Calc
Dynamic Advancing ple ID: culation Model ASTM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse	elp		Calc
Dynamic Advancing ple ID: culation Model ASTM D5946 Polar/Dispersive	SuperLow View Surface Te Select File Water: DIM:	🚰 Browse			Calc




0% Page 1 of 2 PocketGonic

PocketGoniometer ® Program version 2.2

<

1>





Static Dynamic Advancing SuperLow View Surface T	ension Surface Energy Calibrate Help	
Static Dynamic Advancing SuperLow View Surface T	Threshold Scale Result	Start
Start Portectorion et al 100 Porteck DL Microsof The Focketoonnone teles Frugi universion 252 Static Dynamic Advancing SuperLow View Surface T	Threshold Scale Result Tilt Angle:0.2 Degrees Threshold: 140	Next
Insert semi-sphere!		
	Height 0.94 Base 1.99 Volume 1.90 Angle 86.7	